

ESPEC CSH / CWH *Stability Test Chamber*



- suitable for stability test of pharmaceuticals, biologics, goods and personal care products
- meet ICH guideline for stability testing
- support serve storage conditions of $\pm 1^{\circ}\text{C}$, $\pm 5\% \text{ rh}$
- high reliability and easy maintenance
- email alert for faster recovery

Model	Temperature	Humidity	Inside Capacity
CSH (Reach In)	20 to 75 °C	50 to 90 %rh	235L, 470L 794L
CWH (Walk In)	25 to 40 °C	60 to 75 %rh	8 m ³ , 12 m ³ 17 m ³